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OPERATING MANUAL

ELECTRONIC

ABBEREFRACTOMETER AR 2008

INHALT**CONTENTS**

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I. USAGES:

The Abbe-refractometer measures refractive indexes n_D , solid contents in % and mean dispersion values $n_F - n_C$ of transparent and translucent liquids or solids. When connected to a thermostat the instrument can measure the refractive indexes from $0^\circ\text{C} - 70^\circ\text{C}$.

The refractive indexes and the dispersion values belong to the most important optical data of a substance with which the purity, concentration and optical condition can be tested.

The Abbe-refractometer has a wide field of application eg. in the petro-chemistry, oil and lubrication industries and in the pharmaceutical, chemical, food and sugar industries.

It is suitable for use in the following analyses:

- determination of concentration
- determination of mixing ration
- purity control
- operation control of intermediate and end production

II. SETTING UP:

The Abbe-refractometer is delivered in a lockable box. Remove the refractometer and its attachments from the box.

Screw the temperature sensor, with its sealing ring, in the refractometer's thermometer mount (16) and connect the cable to the thermometer socket (19). Screw the lighting fitment, with the knurled screw, onto the lighting mount (22) and attach the connecting cable to the lighting connection (20).

It is necessary to keep the temperature constant for absolute exact measurements because the refractive index of liquids alters with the temperature.

Flexible piping connections with a diameter of 10mm are supplied to keep the prism at a lukewarm temperature.

Connect the thermostat pipe to the connecting nipple (18). The lukewarm water flows round the measuring prism, through a small piece of pipe from nipple (4) to nipple (3), then through the illuminating prism and leaves this through nipple (2).

The instrument is ready for use when the power cable is connected to the power supply socket (6). After switching on the instrument with the mains switch (7) “—> Zero point” always appears on the display screen. After switching on one must always turn the scale adjustment knob (5) clockwise until the zero stop point is

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reached, for the purpose of automatic adjustment of the measuring system. The instrument is then ready for use.

The measuring instrument has an electronic clock showing date and time for the test certificate. Before initial usage of the instrument these functions should be controlled and if necessary corrected.

When the selector button (11) and the confirmation button(10) are pressed simultaneously the instrument can enter the menu. Afterwards the system asks for information which has to be answered with "yes" or "no". To answer the questions move the cursor with the selector button (11) to "Y" for yes or "N" for no. If the question is answered correctly press the confirmation button (10) and the next question will appear on the display.

First of all, the system asks for the type of measurement, the refractive index nD or the solid content in % and if desired with temperature compensation. After answering these questions the system asks "System Y/N". The question must be answered with yes. The instrument can then be programmed to one of four different output records on the serial port, two different languages and asks for date and time.

With the menu "printer" the lookalike of the printing/data transfer of the serial port can be chosen. There are 4 possibilities:

1. print 1 Ger. (for normal DIN-A4-printer):
Messung Nr. 02: 28.01.94 11:12 1,3300nD 20,0°C
Messung Nr. 03: 28.01.94 11:12 temperaturkompensiert 15,4% 20,0°C
Messung Nr. 04: 28.01.94 11:12 15,4% 20,0°C

2. print 2 eng. (for normal DIN-A4-printer):
measure no. 05: 28.01.94 11:15 1,3300nD 20,0°C
measure no. 06: 28.01.94 11:15 temperaturcompensation 15,4% 20,0°C
measure no. 07: 28.01.94 11:15 15,4% 20,0°C

3. print 3 uni. (for data transmission to PC):
28.01.94 11:18 1,3300nD 20,0°C
28.01.94 11:18 T 15,4% 20,0°C
28.01.94 11:18 15,4% 20,0°C

4. print 4 for 24 char printer (always 3 lines):
Measure no.: 8
28.01.94 11:21
1,3300nD 20,0°C (or 15,4% 20,0°C)

If no printout is chosen (All questions answered with "No"), no characters are transmitted to the serial output port!

After answering the questions for the printout the

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language of the display can be chosen. After answering the question for "language" with y the German or English language can be chosen.

Answering the question "time/date" with "Yes" the date and time will be displayed. If this information is incorrect the question "set clock Y/N" will be answered with Yes and the time can then be entered in hours, minutes and seconds and the date in days, months and years. When the selector button (11) is pressed the numbers shown are increased by one, the correct time and date are then displayed when the confirmation button (10) is pressed.

After setting the time and date the instrument returns to the normal measuring mode. Automatic changing from winter to summer time does not take place. With the aid of a small rechargeable battery the electronic clock continues to function when the instrument is not in use.

The chosen parameters are stored permanently, therefore it is not necessary for reprogramming the unit each time it is in use.

III.CALIBRATION:

Before taking any measurements, the calibration of the scale must be checked. Two possible ways of doing this are described as follows :

1. Calibration with standard solids

After switching the instrument on turn the locking knob(15) clockwise and lift up the illuminating prism(23). Put a few drops of bromnaphthalin on the measuring prism before placing the polished surface of the standard glass block on the measuring prism. Make sure that there are no air bubbles between the two surfaces.

The standard glass block must be slightly covered by the illuminating prism(23). This is done so that the direct light is reduced otherwise an over-exposure will occur which would make an exact calibration impossible.

Using the scale adjustment knob(5) set the LCD scale on the refractive index engraved on the calibrating glass block.

Look through the ocular(14) and turn to the left or to the right until it is possible to see the reticle clearly. Then eliminate the colour round the light/dark line with the compensation knob(1), the separation line must now lie on the point of intersection of the reticle. (diagramm 5). If this is not the case, adjust the position of the light/dark line by turning the adjustment screw(13) on the front of the refractometer, with a screwdriver, until it lies exactly on the point of intersection of the reticle.

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
2. Calibration with distilled water

The method described above is theoretically the more exact method but the following method is easier to use and, therefore, is used more frequently.

Instead of using bromnaphtalin, a few drops of distilled water are put onto the measuring prism. Otherwise, the above instructions are to be followed with the exception that a refractive index of 1.3330 nD or 0% Brix is to be set on the LCD scale.

By this method of calibration, it is necessary to take into account that the refractive index of distilled water depends a lot on the temperature; therefore, it is necessary to measure the temperature and with the help of table 1 to set the corrected refractive index on the scale. Then, as above, the light/dark line can be set on the point of intersection of the reticle.

3. Calibration of the temperature display

The unit is equipped with an electronic temperature sensor. After some time or after shocking the unit it may be necessary for recalibrating the temperature display. Because the temperature is part of the temperature-compensation, a wrong temperature may alter the result of the measurement. For calibrating the temperature a second precision thermometer is necessary.

Before calibrating the ground plate of the instrument is to be removed. There are three adjustment parts (see next illustration) on the electronic board:

The temperature probe has to be dismantled from the unit. The probe should be put in ice-water and with the knob 1 the display is calibrated to 0.5°C. Afterwards the probe together with the probe of the second thermometer is putted in warm water. With knob 2 the display is calibrated to the display of the second thermometer. In both cases the probes should for 15 minutes in the water before calibrating. Now the instrument is calibrated and the ground plate can be mounted to the unit. Knob 3 is for selecting the contrast of the LCD-display.

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IV. TO MEASURE LIQUIDS:

Initially, it must be stipulated whether the results should be shown as refractive index nD values or as solid contents in % Brix values. First of all, push the two buttons (10) and (11) simultaneously. The display shows “Refractive index nD Y/N”, whereby the cursor blinks under the “N”. The cursor can be moved between “Y” and “N” by means of the selector button (11).

If one answers with “N” and pushes the confirmation button (10) the next question that appears is “Brix % Y/N”. When answered with “Y” by moving the cursor and pressing the confirmation button, the system then asks if the result should be shown with or without temperature compensation. After answering these questions the question “System Y/N” will be shown on the display. This question will normally be answered with “N” and the instrument then returns to the measuring mode.

Turn the locking knob (15) to lift up the illuminating prism (23).

Put one or two drops of liquid on the measuring prism, then lower the illuminating prism to its former position and lock with the knob (15).

Look through the ocular (14) and adjust the brightness with the brightness control knob (17) and also sharply define the reticle by turning to the right or left.

By turning the scale adjustment knob (5) move along the measuring range until a light/dark division appears. Get a sharp division line with the compensation knob (1) until a black-white division is reached.

Adjust the light/dark line exactly on the point of intersection of the reticle and then the display will show the refractive index or the solid content in %.

At the same time, the display shows the present temperature of the measuring prism. By pressing the confirmation button (10), the measurement values together with the date and time can be printed on a printer that is connected to the data output (9). If a computer is connected to the data output (8) or (9) then the measurement values plus date and time will be passed to it when the confirmation button (10) is pressed.

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V. MEASUREMENT OF SOLID SUBSTANCES:

The surface of the object to be measured should have a plane smooth surface.

Put one or two drops of bromnaphthalin on the measuring prism and lightly press the smooth surface of the sample onto it. There must be no air bubbles. Open the reflecting mirror (21) and direct light by way of the mirror onto the prism. Dim this light by means of the brightness control knob (17).

Proceed in the same way as for measuring liquids.

The measurements occur now in reflected light, the light/dark fields are exchanged and the contrast becomes less but this has no effect on the measurement results.

When measuring solids the refractive index of the contact fluid must be higher than that of the solid to be measured. When the refractive index of the solid to be measured is higher than 1.65 nD then methyljod with 1.71 nD must be used as a contact fluid.

The surfaces of the prisms must be cleaned immediately after use with alcohol because bromnaphthalin and methyljod are very corrosive.

The prisms are made of reasonably soft material, therefore, one must be careful not to scratch the surfaces.

VI. DETERMINATION OF DISPERSION VALUES:

The dispersion value can be calculated from the following equation:-

$$n_F - n_C = A + B + a$$

The values for A & B in table 2 are related to refractive indices. The values in between are to be obtained by interpolation.

The refractive index of the sample must be measured, the values appertaining to A & B are to be taken from table 2.

The alpha value is determined as follows :

When measuring the refractive index, the colour on the dividing line is eliminated by using the compensation knob (1). There is a scale ring on the knob which is divided into 60-0-60.

There are two positions on the compensation knob where it is possible to eliminate the colour fringe on the dividing line.

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Adjust to the two positions and read the corresponding value on the scale ring on the compensation knob(1). This process should be repeated several times and the values noted on a table so that the mean value can be calculated.

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The alpha value can then be determined using the mean value “Z” according to table 2. The values in between can be found by interpolation.

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The Abbe figure is another important quantity as an identification of an optical medium. It is calculated as follows:

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Calculation:

n _D	=	1.4827	(measured)
Z	=	42.2	(mean)
A	=	0.024265	(from table)
B	=	0.024670	(from table)
a	=	-0.5962	(from table)
B x a	=	-0.014708	
n _F - n _C	=	A + B x a	
	=	0.009557	

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VII. SPECIAL INSTRUCTIONS:

The contrast of the light/dark area is less if the sample is cloudy or coloured. In this case the measurement must be done with reflected light as described for solid samples.

When a colloid solution is used one must make sure that the colloids are finely and homogeneously distributed in the solution otherwise the dividing line will appear to be blurred.

A highly viscous solution could show a different refractive index on the surface as on the bottom. The solutions should always be well stirred and mixed.

The solution must cover the larger part of the measuring prism otherwise no clear light/dark division can be obtained. This can, in particular, happen by quickly evaporating liquids that remain a longer time between the measuring and the illuminating prisms.

To have optimal adhesion of the samples the prism must always be kept clean.

The instrument must not come in contact with dust, gases, direct sunlight or dampness. It must always be used in an horizontal position otherwise incorrect values will be measured.

VIII. REPLACEMENT OF THE MEASURING PRISM:

After extensive use the polished surface of the measuring prism often becomes corroded.

Therefore, it must be replaced by a new one. The fitment with the measuring prism can easily be loosened and removed.

The thermostat tube (4) must be unscrewed counter-clockwise and the nut of the prism support removed with the special spanner. The measuring prism can then be exchanged. Afterwards the instrument must be readjusted.

The order number for the replacement with the special spanner is AR 11.

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IX. TECHNICAL DATA:

measuring range :
 refracting index nD 1.3000-1.7200
 solid content in % 0.0-95.0

resolution and accuracy :
 refracting index nD 0.0001/0.0002
 solid content in % 0.1/0.2

main supply : 115/230 V AC, 50 Hz, 50 W

lightning : light emitting diode with 589 nm

measurement of temperature : digital 0° to 99° C




dimensions : 140x100x235 mm

weight : 5 kg / 11 lbs

CE: The instrument corresponds to these rules:
 EN 50081-1, IEC 1000-4-2, IEC 1000-4-4, IEC 1000-4-5

dataports :
 connector 9 seriell RS 232 9600 baud, 8 bit,
 1 stopbit, no parity
 connector 8 seriell RS 422 9600 baud, 8 bit,
 1 stopbit, no parity

Connections :
 port 9 (RS 232): TX-data on 3, RX-data on 2,
 (Diagram 3) ground on 7
 port 8 (RS 422): TX-data + on 2, TX-data - on 9,
 (Diagram 4) ground on 8

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
X.TWOPPOINT CALIBRATION

Once a year or when misreadings are recognized, the main calibration of the instrument should be checked. This main calibration is a two - point - calibration, this means, the instrument is calibrated to a lower and to an upper reading.

With the AR 2008 this two calibrations are done for 0 % Brix = 1.3330 nD R.I. and the upper point for 1.5061 nD R.I.. The lower point is calibrated with distilled water (must be done at 20°C regarding to the temperature compensation) and the upper point is calibrated with the standard glass block supplied with the instrument.

Testing this two calibration points and getting proper readings, no other work has to be done. Otherwise a new two-point-calibration has to be done by switching the instrument to the service mode. Switch the instrument off and remove the lower plate.

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On the CPU board put DIP-switch no 8 to "ON" position. Afterwards switch instrument on and turn the scale adjustment knob (5) clockwise until the stop point is reached.

Put some drops distilled water on the measuring prism and close and lock the illumination prism(23). Look through the ocular(14) and adjust the light/dark line to the recticle with the scale adjustment knob. Push the right key(10). The lower calibration point is set and stored.

Open the illumination prism, remove the water and position the standard glass block with bromnaphthalin on the measuring prism (see chapter calibration with standard solids). Look through the ocular(14) and adjust the light/dark line to the recticle with the scale adjustment knob. Push the right key(10). The upper calibration point is set and stored.

Switch off the instrument. switch back DIP-switch 8 to "OFF" and fix the lower plate. The instrument is calibrated.

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
XI. FLOWING THROUGH CELL

The instrument can be equipped with a flowing through cell. There is a flowing through cell with a funnel AR 15 or a flowing through cell for a pump AR 16 available. The sample passes the prism continuously with both flowing through cells. With the flowing through cell AR 16 a pump must be used.

The procedure for connecting the flowing through cell is the same for both models.

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Probe Nr. : 10
26.05.98 17:08
TC 7,3% 21,3°C

XII. PRINTER

You can connect the instrument to the **printer CBM 910**. The readings will be printed with date, time and sample number.

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You connect the printer with the included cable to the serial port (9) of the instrument. In the system menu in menu printer the function print 4 must be chosen (see page 5). Pushing **PRINT key (10)** prints the reading with date, time and sample number. After switching on, the instruments counts the sample number from 1 upwards by itself.

Table 1 :

Refracting index and dispersion of distilled water by means of temperature.

°C	nD	nF - nC	°C	nD	nF - nC
10	1,33369	0,00600	26	1,33240	0,00596
11	1,33364	0,00600	27	1,33229	0,00595
12	1,33358	0,00599	28	1,33217	0,00595
13	1,33352	0,00599	29	1,33206	0,00594
14	1,33346	0,00599	30	1,33194	0,00594
15	1,33339	0,00599	31	1,33182	0,00594
16	1,33331	0,00598	32	1,33170	0,00593
17	1,33324	0,00598	33	1,33157	0,00593
18	1,33316	0,00598	34	1,33144	0,00593
19	1,33307	0,00597	35	1,33131	0,00592
20	1,33299	0,00597	36	1,33117	0,00592
21	1,33290	0,00597	37	1,33104	0,00591
22	1,33280	0,00597	38	1,33090	0,00591
23	1,33271	0,00596	39	1,33075	0,00591
24	1,33261	0,00596	40	1,33061	0,00590
25	1,33250	0,00596			

Table 2 :

Table of dispersion for ABBE-refractometer :

Calculation: $n_D - n_C = A + B \times a$

For Z = 0 .. 30 plus (+) for a

For Z = 30 .. 60 minus (-) for a

nD	A	B	Z	a	Z
1,30000	0,02496	0,02901	0	0,000	60
1,31000	0,02490	0,02889	1	0,999	59
1,32000	0,02485	0,02876	2	0,995	58
1,33000	0,02480	0,02861	3	0,988	57
1,34000	0,02476	0,02845	4	0,978	56
1,35000	0,02471	0,02828	5	0,966	55
1,36000	0,02466	0,02810	6	0,951	54
1,37000	0,02462	0,02790	7	0,934	53
1,38000	0,02458	0,02768	8	0,914	52
1,39000	0,02454	0,02746	9	0,891	51
1,40000	0,02450	0,02722	10	0,866	50
1,41000	0,02446	0,02696	11	0,839	49
1,42000	0,02443	0,02670	12	0,809	48
1,43000	0,02440	0,02641	13	0,777	47
1,44000	0,02437	0,02612	14	0,743	46
1,45000	0,02434	0,02580	15	0,707	45
1,46000	0,02431	0,02547	16	0,669	44
1,47000	0,02429	0,02513	17	0,629	43
1,48000	0,02427	0,02477	18	0,588	42
1,49000	0,02425	0,02440	19	0,545	41
1,50000	0,02423	0,02400	20	0,500	40
1,51000	0,02422	0,02359	21	0,454	39
1,52000	0,02421	0,02316	22	0,407	38
1,53000	0,02420	0,02272	23	0,358	37
1,54000	0,02419	0,02225	24	0,309	36
1,55000	0,02419	0,02176	25	0,259	35
1,56000	0,02419	0,02126	26	0,208	34
1,57000	0,02420	0,02073	27	0,156	33
1,58000	0,02421	0,02017	28	0,104	32
1,59000	0,02422	0,01959	29	0,052	31
1,60000	0,02424	0,01899	30	0,000	30
1,61000	0,02427	0,01835			
1,62000	0,02430	0,01769			
1,63000	0,02434	0,01699			
1,64000	0,02438	0,01626			
1,65000	0,02443	0,01549			
1,66000	0,02450	0,01467			
1,67000	0,02457	0,01381			
1,68000	0,02466	0,01288			
1,69000	0,02476	0,01189			
1,70000	0,02489	0,01083			