

MS-BV 1-1000

Stainless steel discs measuring system
for viscosity measurements in a 150 ml Beaker.
This system allows rapid and economical measurements.
This system could be used by First RM
and can also be used with the Rheomat RM 100, RM 200.

measuring system

	Name	Reference	Dimensions	Volume	Viscosity range
	BV 1-100 axis	117102	-	-	-
	BV centring device	117202	Ø 46.93 mm	-	-
	BV disc n°1	117001	Ø 45 mm	120 ml	2 - 500 000 mPa.s
	BV disc n°10	117010	Ø 40 mm	120 ml	20 - 5 000 000 mPa.s
	BV disc n°100	117100	Ø 20 mm	120 ml	200 - 45 000 000 mPa.s
	BV 1000 axis	117101	Ø 4 mm	120 ml	2000 - 510 000 000 mPa.s
	Glass Beaker 150 ml	117150	Ø 50-52 mm	150 ml	-
	MS TI tube	118001	Ø 50 mm	150 ml	-

APPLICATIONS

Food industry



Cosmetics/Pharmaceuticals



Paint



Chemical industry

