

MS-BV 1-1000

[← PRODUCTS LIST](#)

Stainless steel discs measuring system for viscosity measurements in a 150 ml Beaker.

This system allows rapid and economical measurements.

This system could be used by First RM and can also be used with the Rheomat RM 100, RM 200.

[REQUEST A QUOTE](#)[DOWNLOAD THE BROCHURE](#) [ZOOM](#)

Name	Reference	Dimensions	Volume	Viscosity range
BV 1-100 axis	117102	-	-	-
BV centring device	117202	Ø 46.93 mm	-	-
BV disc n°1	117001	Ø 45 mm	120 ml	2 - 500 000 mPa.s
BV disc n°10	117010	Ø 40 mm	120 ml	20 - 5 000 000 mPa.s
BV disc n°100	117100	Ø 20 mm	120 ml	200 - 45 000 000 mPa.s
BV 1000 axis	117101	Ø 4 mm	120 ml	2000 - 510 000 000 mPa.s
Glass Beaker 150 ml	117150	Ø 50 - 52 mm	150 ml	-
Polypro. Becher 150 ml	117151	Ø 56 mm	150 ml	-
MS TI tube	118001	Ø 50 mm	150 ml	-